

CATS Component and Activity Tracking System

A substantial amount of labor and downtime is associated with routine maintenance, repair, and tracking of materials and consumables. Gas cylinder management and tracking recertifications is time-consuming as well, and errors in managing all of these components can cast previously monitored data into question.

And then there's all the paper forms- there's got to be a better way!

Using Agilaire's new **Component and Activity Tracking System (CATS)** and a wireless device (or even your mobile phone), each field technician can remotely access your parts and consumable inventory, locate items, mark their location and usage, and maintain your gas cylinder inventory through easy-to-use forms. Synchronize instantly from the central database to keep everything and everybody up to date. Use optional barcode labels for advanced tracking of bottles and parts usage.

Task	
Assigned to User:	Will Peters
Type:	QC_Check
Work Location:	
Date Due:	<input checked="" type="checkbox"/> 2/ 4/2009
Repeat Int Days: ?	91
Date Completed:	<input type="checkbox"/> 2/25/2009
Description:	
Quarterly Audit for Ozone	
<div style="border: 1px solid black; height: 100px; width: 100%;"></div>	
<div style="display: flex; justify-content: space-between;"> Task Info Consumables Parts </div>	

Easily create historical logs of activity and usage, and monitor for low inventory conditions.

Track long-term costs for each analyte.



Your day is stretched as it is.

Agilaire helps you breathe easier with CATS / CATS-PAW - a software solution for site activity and data management.

Finally, you can breathe easier about your inventory!

For more information, logon at: www.agilairecorp.com/cats

CATS Component and Activity Tracking System

Client Requirements

- Mobile Windows Device
- Optional: Barcode Scanner, Digital Camera
- Wireless data plan

Server Requirements


- Windows IIS (Internet Information Server) with .NET 3.5
- SQL Express or Workgroups

Capabilities

- Instrument / Part Templates
- Individual Component Tracking
 - Denote Analyzer Condition, Retirement
- Track Consumable Usage
- Track Instrument Repairs
- Track Bottle Expirations and Recertification Events
- Track Bottle Pressure Inspections
- Activity Tracking by User, Due Date, Site
 - Designate Parts/Consumables For Activities
- Schedule Repeating Tasks with Email Notification
- Designate Inventory Sites
- Server-based Summary Reports
- Purchase Order Tracking / Optional Purchasing System Integration
 - "On Order" vs. "Received" Status

Features

- Bidirectional Data Synchronization
- Input Forms with pull-down pick lists
- Signature Blocks
- Camera Annotation
- Comment Fields with Each Record



For expanded capability, use our Portable Audit Wizard for CATS (CATS-PAW) and eliminate the paper headache (and transcription errors) from your field audits. Agilaire offers customization of templates to meet individual agency requirements.